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Rx Daytona Plus:

System Setup 1

Host Communication Mode
 Off Line
 On Line Batch
 On Line Real time

Host Details Mode
 Auto rerun
 Confirm host for rerun
 Inquire host for rerun

Host Communications
Baud Rate: 9600
Data Bit: 8bit
Stop Bit: 1bit
Parity Bit: Even
Protocol (E13B1): Normal

Others
 Disable extrapolated results
 Auto scan for Reagent Container
 Check remaining reagent volume
 Lamp voltage control off
 Initial Result of Rerun
 Auto Rerun
 ~~Print test~~ sequence number
 Enable Auto Backup
 Book orders for future rounds

Live Print out
 Enable
Report
 Enable Live QC Report

Printout Header

1
2
3
4
5


Prepričajte se, da to polje ni označeno.

Rx Imola:

The screenshot displays the 'Run Monitor' software interface. At the top, there are tabs for 'Run (F5)', 'Parameter (F6)', 'Calibration (F7)', 'QC (F8)', 'System (F9)', 'Mainte (F10)', 'EL', and 'SHUTDOWN'. The main area is divided into a left control panel and a central circular well plate visualization. The control panel includes fields for 'Time at run started', 'Estimated time at sampling complete', and 'Estimated time at run end'. It also features a 'Halogen' setting at 20 Min and an 'Auto Rerun' checkbox, which is circled in red. Below these are buttons for 'SPT Wash', 'MIX Ascend', 'WU Ascend', 'Inner ASP', 'Initialization', 'Prime' (set to 1), 'ISE Prime' (set to 5), 'ISE Calibration', and 'ISE Cleaning'. The central well plate shows 72 numbered wells, each with a status indicator. A legend on the left explains the status colors: Green (Sampling Started), Blue (Range Over), Purple (Rerun Required), Red (Error), Cyan (Process Complete), White (Not Processed), and Yellow (No test Ordered). The 'Indication' field at the bottom shows 'Automatic Rerun.' and there are buttons for 'Start (F1)', 'Stop (F2)', 'STAT (F3)', 'Alarm (F4)', 'TabUP (F11)', and 'TabDOWN (F12)'.

Prepričajte se, da to polje ni označeno.

Rx Modena:

Prepričajte se, da to polje ni označeno.