VARNOSTNO OBVESTILO

DiaClon Anti-P1, Ref. 007212

Pomembne informacije za DiaClon Anti-P1, Ref. 007212

V vednost profesionalnim uporabnikom v laboratoriju.

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B007212 v.02 2024-08-03 / Napačna različica		
Well Reaction Grade	Reaction Description	Result Interpretation
-	A flat and compact pellet of RBCs with a smooth surface at the bottom of the well with no visible agglutination.	Negative (indicates an absence of the corresponding antigen)
+/-	A pellet of RBCs at the bottom of the well with very few agglutinated RBCs visible above the pellet or an irregular pellet.	Non-interpretable (further investigations are required)
+	A pellet of RBCs at the bottom of the well with agglutinated RBCs visible in the lower half of the gel column.	
++	Agglutinated RBCs distributed throughout the entire length of the gel column, with no line of RBCs on the top of the well.	
+++	Most agglutinated RBCs concentrated at the top of the gel or upper half of the gel column.	Positive (indicates the presence of the corresponding antigen)
++++	Agglutinated RBCs concentrated as a line on the top of the gel column with a few agglutinated RBCs just underneath the gel surface.	
dp: double population	Agglutinated RBCs as a line at the top of the gel or dispersed in the upper part of the gel and non-agglutinated RBCs forming a pellet at the bottom of the well.	Non-interpretable (further investigations are required)

B007212 v.02.1 2025-11-04 / Popravljen različica		
Well Reaction Grade	Reaction Description	Result Interpretation
-	A flat and compact pellet of RBCs with a smooth surface at the bottom of the well with no visible agglutination.	Negative (indicates an absence of the corresponding antigen)
+/-	A pellet of RBCs at the bottom of the well with very few agglutinated RBCs visible above the pellet or an irregular pellet.	Non-interpretable (further investigations are required)
+	A pellet of RBCs at the bottom of the well with agglutinated RBCs visible in the lower half of the gel column.	Positive (indicates the presence of the corresponding antigen)
++	Agglutinated RBCs distributed throughout the entire length of the gel column, with no line of RBCs on the top of the well.	
+++	Most agglutinated RBCs concentrated at the top of the gel or upper half of the gel column.	
++++	Agglutinated RBCs concentrated as a line on the top of the gel column with a few agglutinated RBCs just underneath the gel surface.	
dp: double population	Agglutinated RBCs as a line at the top of the gel or dispersed in the upper part of the gel and non-agglutinated RBCs forming a pellet at the bottom of the well.	Non-interpretable (further investigations are required)

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Elizabeth Platt

VP, Quality, Regulatory & Clinical Affairs
Bio-Rad Laboratories Inc.

Enclosure:
FIELD ACTION RESPONSE FORM

FIELD ACTION RESPONSE FORM
Bio-Rad Reference: FA 2025-05
Bio-Rad Product Segment: IH
Single Registration Number (SRN): CH-MF-000020826

PRODUCT

DiaClon Anti-P1, Catalog Reference 007212
<i>e-IFU available since October 09th, 2025 in downloads.bio-rad.com: e-IFU code number B007212 v.02 2024-08-03</i>

CUSTOMER / CHANNEL PARTNER INFORMATION

Account Name:	
Undersigning Manager Name:	
Address:	
Telephone Number / Fax:	
Account Number:	

STATEMENT:

- No affected product received
- I am aware of the information about the field action concerning the above reference product(s) and have proceeded according to the instructions issued by Bio-Rad.
- <Omit/include, as needed>* For completion by Channel Partners: All customers have been informed about this field action and have proceeded according to the instructions issued by Bio-Rad. Number of customers informed: _____

Number of affected products received:	N/A	Number of affected products corrected/ destroyed/ returned (as applicable to the Field Action instructions):	N/A
If number of products corrected/ destroyed/ returned is different to the number received, please account for the difference: N/A			

Date:
 Customer / Channel Partner Signature (and Stamp if applicable)
 Please return this form to: **info@diahem.si**