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**DiaClon Anti-C<sup>W</sup>** 4/7  
English B001321 v.02 2025-04-29

**VISUAL READING AND INTERPRETATION OF THE RESULTS**

Well Reaction Grade	Reaction Description	Result Interpretation
-	A flat and compact pellet of RBCs with a smooth surface at the bottom of the well with no visible agglutination.	Negative (indicates an absence of the corresponding antigen)
+/-	A pellet of RBCs at the bottom of the well with very few agglutinated RBCs visible above the pellet or an irregular pellet.	Non-interpretable (further investigations are required)
+	A pellet of RBCs at the bottom of the well with agglutinated RBCs visible in the lower half of the gel column.	Non-interpretable (further investigations are required)
++	Agglutinated RBCs distributed throughout the entire length of the gel column, with no line of RBCs on the top of the well.	Positive (indicates the presence of the corresponding antigen)
+++	Most agglutinated RBCs concentrated at the top of the gel or upper half of the gel column.	Positive (indicates the presence of the corresponding antigen)
++++	Agglutinated RBCs concentrated as a line on the top of the gel column with a few agglutinated RBCs just underneath the gel surface.	Positive (indicates the presence of the corresponding antigen)
dp: double population	Agglutinated RBCs as a line at the top of the gel or dispersed in the upper part of the gel and non-agglutinated RBCs forming a pellet at the bottom of the well.	Non-interpretable (further investigations are required)

Non-interpretable results on instruments require edition by an operator.  
Any other reactions not in accordance with the table above are non-interpretable and require further investigations.  
In case a positive reaction is obtained, the test must be repeated, and the ID-Card "Control Card A" (Id-n° 52020) processed in parallel according to the "TEST PROCEDURE".  
It is recommended to include the "Control Card A" (Id-n° 52020) when using the ID-Card "DiaClon Anti-C<sup>W</sup>". A negative reaction obtained with the "Control Card A" validates the results obtained by using the "DiaClon Anti-C<sup>W</sup>". In case a positive reaction is obtained with the "Control Card A" (Id-n° 52020), the antigen determination is not valid. Further investigation is required.

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Bio-Rad and the Bio-Rad logo are trademarks of Bio-Rad Laboratories, Inc. in certain jurisdictions.