Medline International Germany GmbH – Medline Str. 1-3 – D-47533 Kleve

Company Name

Address

Address

ZIP City

Kleve, 29. Juli 2026

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




REF ORNEX610VCA


en MED-STOP Vacuum Control Valve

OR853VCAMED-STOP, Adult Manual Vacuum Controller
 OR854VCPMED-STOP, Pediatric Manual Vacuum Controller
 ORNEX66VCAORNEX Suction Tubing (1.8m) + Adult MED-STOP Control Valve
 ORNEX66VCPORNEX Suction Tubing (1.8m) + Pediatric MED-STOP Control Valve
 ORNEX610VCAORNEX Suction Tubing (3m) + Adult MED-STOP Control Valve

Indication for Use:

1. MED-STOP Vacuum Control Valve Setup


- Connect the MED-STOP Vacuum Control Valve to the suction tubing.
- Proceed to rinse the device by suctioning sterile saline or sterile water to verify device functionality to ensure proper operation before patient application.
- MED-STOP is ready for use.

2. MED-STOP Vacuum Control Valve Use

- Connect the MED-STOP Vacuum Control Valve to the suction catheter and proceed with suction.
- To manually apply suction, cover the open hole on the blue cover with a finger.
- When the suction system and connected MED-STOP are not in use during a procedure, ensure the vacuum is turned off.
- Refer to rinsing during use section to prevent blockage from suctioned material.

3. MED-STOP Vacuum Control Valve Rinsing During Use

- Carefully remove the suction catheter and insert the MED-STOP Vacuum Control Valve in the cleaning bottle cap.
- Rinse the device during use by suctioning sterile saline or sterile water, as clinically indicated, to maintain system patency and prevent blockage from secretions, mucus, or other material. Replace the device if obstructions cannot be cleared.

No need to follow STEP 1 for ORNEX + MED-STOP set (ORNEX66VCA, ORNEX66VCP and ORNEX610VCA).

Contraindications:

Do NOT use if the valve does not work correctly.

Important:

- Reuse of single use devices creates a potential risk for the patient or user. It may cause contamination and/or impairment of function, which may lead to injury, illness and/or death.
- If a serious incident occurs in relation to the use of the device, it should be reported to Medline and relevant competent authorities of the Member State in which the user and/or patient is established.
- Single-use device — Do not reuse. Discard the device after the procedure in accordance with facility protocol and use a new MED-STOP device for each new procedure.
- Intended for use during single suction procedures.
- Use only by trained healthcare professionals familiar with suction equipment.
- Prior to use, visually inspect valve for obstructions or particulate. Do not use the product if obstructions or particulate matter are present; discard the product in accordance with local procedures.
- Prior to use, rinse the device by suctioning sterile saline or sterile water to verify device functionality to ensure proper operation before patient application.
- Rinse the device during use by suctioning sterile saline or sterile water, as clinically indicated, to maintain system patency and prevent blockage from secretions, mucus, or other material. Replace the device if obstructions cannot be cleared.
- When the suction system and connected MED-STOP are not in use during a procedure, ensure the vacuum is turned off.

ORNEX + MED-STOP (ORNEX66VCA, ORNEX66VCP, ORNEX610VCA)

DEHP warning: This product contains the phthalate DEHP. Phthalates are widely used as plasticizers in medical products and have undergone extensive testing. However, when using the device, it should be taken into consideration that possible health effects when used for long durations in pregnant women, nursing mothers and children cannot be entirely discounted.

Instructions for Disposal:

Take precautions when discarding this device. The device should be disposed in accordance with applicable hospital or national regulations for biologically hazardous waste.

Intended Purpose:

Medline Suction Tubing Family products are intended to be used in various medical and surgical procedures to sample, aspirate or to remove body fluid.

Intended Users :

Intended to be used by healthcare professionals and are not intended to be used outside of the supervision of a clinically trained professional. Intended to be used in healthcare environment, including acute care, long-term acute care, post-acute care, long term care/nursing home, surgical/operating room, and procedural room.

Intended Patient Population:

Intended for use on pediatric and adult patient populations, undergoing a procedure to sample, aspirate, or remove bodily fluids.

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Medline International France SAS,
Medline Industries Ltd,
Medline International B.V.,