«Hospital_Name»
«Users_Name»
«Department»
«Customer_Address»
«Zip_Code» «City»
«Country_name»

<Referenca: 97518004-FA>

SRN: IL-MF-000016033

30. marec 2026

Nujno varnostno obvestilo – obvestilo o izdelku Servisne dejavnosti za CO₂ laserske sisteme in sisteme z dvojno valovno dolžino (DW) VersaPulse PowerSuite (VPPS)

Spoštovani «Users_Name»,

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Dodatek 1 Revidirana navodila za servisiranje VPPS

Povzetek sprememb: posodobite v servisnega priročnika VPPS (0636-499-00)	Posodobitev servisnega priročnika: Razdelek 3. Kalibracija, nastavitve in poravnava: 3.0.1 Pregled	Dodatne informacije																		
	POTREBNO ORODJE IN OPREMA: orodje za poravnavo tarče SP0638-497-51	Dodajte orodje za poravnavo tarče kot zahtevano opremo																		
	<p>3.2.2.3 Poravnava resonatorja: koraki 4(h)–4(m)</p> <p>h. Remove the Cross-Hair Target tool and insert the Target Alignment tool (P/N SP0638-497-51) for adjusting the OC.</p> <p>i. Place a chemical Burn Paper (3207-0091) for testing the laser footprint and make a burn spot on the paper. If the mark is not solid and symmetrical, adjust the OC adjustment screws to achieve a nice, round pattern per the example below:</p> <div style="text-align: center;"> <table border="1" style="margin: auto;"> <tr> <td style="padding: 5px;">Good Sample</td> <td style="padding: 5px;">Bad Sample</td> </tr> <tr> <td style="text-align: center;"></td> <td style="text-align: center;"></td> </tr> </table> </div> <p>j. Repeat steps e-i until the spot size is centered in the Target Alignment tool and a clear, round, symmetric spot is achieved on the burn paper.</p> <p>k. After aligning the resonator, tighten the screws and check that the spot is still centered in the Target Alignment tool, and a clear, round, symmetric spot is achieved on the burn paper.</p> <p>l. Once the burn mark is centered within the crosshairs and the spot size is centered in the Target Alignment tool, remove the tool and position the power meter head behind the aperture so that the YAG energy passing through the hole is directed onto the power meter sensor. Activate the laser and measure the output power in accordance with the table below. Adjust the lamp energy/current as needed to reach the specified minimum power (for Nd:YAG, begin at 10 A).</p> <table border="1" style="margin: auto; width: 100%; border-collapse: collapse;"> <thead> <tr> <th style="text-align: left;">YAG</th> <th style="text-align: left;">Lamp voltage</th> <th style="text-align: left;">Lamp Max energy/current</th> <th style="text-align: left;">Frequency</th> <th style="text-align: left;">Required minimum power</th> </tr> </thead> <tbody> <tr> <td>HO:YAG</td> <td style="text-align: center;">800v</td> <td style="text-align: center;">75J</td> <td style="text-align: center;">12.5 Hz</td> <td style="text-align: center;">34W</td> </tr> <tr> <td>Nd:YAG</td> <td style="text-align: center;">800v</td> <td style="text-align: center;">60A</td> <td style="text-align: center;">60 Hz</td> <td style="text-align: center;">70W</td> </tr> </tbody> </table> <p>m. Reinstall the first Relay Mirror and perform the Fiber Alignment Procedure in Topic 3.1.4.</p>	Good Sample	Bad Sample			YAG	Lamp voltage	Lamp Max energy/current	Frequency	Required minimum power	HO:YAG	800v	75J	12.5 Hz	34W	Nd:YAG	800v	60A	60 Hz	70W
Good Sample	Bad Sample																			
																				
YAG	Lamp voltage	Lamp Max energy/current	Frequency	Required minimum power																
HO:YAG	800v	75J	12.5 Hz	34W																
Nd:YAG	800v	60A	60 Hz	70W																

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