



April 10, 2026

Ref #: FA-2026-020 EMEA

URGENT: VOLUNTARY PRODUCT FIELD ACTION

RE: O&M Halyard, Inc. Product Field Action – HALYARD* PUREZERO* Ultraviolet* Nitrile Exam Gloves

Dear Valued Customer,

A Product Field Action is being initiated by O&M Halyard, Inc. related to the HALYARD* PUREZERO* Ultraviolet* Nitrile Exam Gloves (codes 48771 -48775). This Field Action is related to not meeting the minimum breakthrough time of 30 minutes for n-Heptane required for European Type B performance per ISO 374-1:2016+A1:2018 and EN 16523-1:2015+A1:2018. While O&M Halyard, Inc. has received zero complaints or reports of adverse events that would indicate a real-world concern, O&M Halyard, Inc. is issuing this voluntary Field Action.

EVENT DETAILS:

This action affects the product codes identified in the table below. Table 1A identifies the lots distributed between September 2025 and March 2026. O&M Halyard, Inc. records indicate you have received HALYARD* PUREZERO* Ultraviolet* Nitrile Exam Gloves from one of these specific lots that is associated with this action as listed below. The issue surrounding this Field Action is related to not meeting the minimum breakthrough time of 30 minutes for n-Heptane required for European Type B performance per ISO 374-1:2016+A1:2018 and EN 16523-1:2015+A1:2018. No other product codes or lots of gloves are impacted in this Field Action.

Affected Products:

| Product Name | Alternative Code | Code |
|---|------------------|-------|
| HALYARD* PUREZERO* Ultraviolet* Nitrile Exam Gloves | LFS511XS | 48771 |
| HALYARD* PUREZERO* Ultraviolet* Nitrile Exam Gloves | LFS511SM | 48772 |
| HALYARD* PUREZERO* Ultraviolet* Nitrile Exam Gloves | LFS511MD | 48773 |
| HALYARD* PUREZERO* Ultraviolet* Nitrile Exam Gloves | LFS511LG | 48774 |
| HALYARD* PUREZERO* Ultraviolet* Nitrile Exam Gloves | LFS511XL | 48775 |

Refer to the attached table 1A for the specific lots impacted.

Risk to Health: The associated risk with the use of these exam gloves is potential skin/tissue damage or irritation.

Actions to Take:

Please review your inventory, segregate and quarantine all affected product(s). All other product codes and lots of this product are NOT impacted by this Field Action.

Distributors:

- REVIEW** your inventory for the affected **product code(s)** and **specific lot number(s)**.
- Contact the Customer Service team and send email to OMRA_RECALLS@owens-minor.com to request a return goods authorization (RGA) for the affected products.
- COMMUNICATE** this voluntary Field Action to your staff to ensure awareness of it.
- If any of the affected lot(s) identified above were further distributed, please ensure all end-users are appropriately notified of this Field Action and maintain records of effectiveness.
- Any adverse events associated with the use of the affected product should be reported to complaints@owens-minor.com.



6. RETURN the completed Field Action Response Form to OMRA_RECALLS@owens-minor.com as soon as possible even if you have no current product inventory.
7. If credit is requested, please contact the Customer Service team and send email to GM-OMRA_RECALLS@owens-minor.com.

If you have any questions, please email OMRA_RECALLS@owens-minor.com using the Event # FA-2026-020 in the header of the response. This will ensure prompt response by the O&M Halyard, Inc, team.

Thank you,

Caitlin Senter, MS, RAC
Director, Global Regulatory Affairs
O&M Halyard, Inc.

Enclosures:
Field Action Response Form



Table 1A - Specific Lots Impacted

| Product Name | Alternate Code | Code | Specific Lots Impacted |
|---|----------------|-------|------------------------|
| HALYARD* PUREZERO* Ultraviolet* Nitrile Exam Gloves | LFS511XS | 48771 | SJ408532X |
| HALYARD* PUREZERO* Ultraviolet* Nitrile Exam Gloves | LFS511SM | 48772 | SJ415532X |
| | | | SJ428132X |
| HALYARD* PUREZERO* Ultraviolet* Nitrile Exam Gloves | LFS511MD | 48773 | SJ421132X |
| | | | SJ430202X |
| | | | SJ431604X |
| HALYARD* PUREZERO* Ultraviolet* Nitrile Exam Gloves | LFS511LG | 48774 | SJ413433X |
| | | | SJ428832X |
| | | | SJ432304X |
| HALYARD* PUREZERO* Ultraviolet* Nitrile Exam Gloves | LFS511XL | 48775 | SJ413432X |